


Liver vein segmentation for contrast-enhanced MR images

László Ruskó, GE Healthcare, Szeged

Imagination at work.

Overview


- Motivations
- Liver segmentation
- Vessel segmentation
- Challenges of vessel segmentation
- Improvement opportunities




Motivation

Ultrasound – MRI image registration


US-guided radiofrequency ablation




Landmarks: liver vessels





Pre-operative MRI





Intra-operative US




Fusion

Examples of automated registration



Liver segmentation



Liver segmentation



Characteristics of EOB-enhanced MR LAVA

Contrast-enhanced CT





Liver (red) intensity is represented by the last peak of the image histogram (blue), and it is always found in a well defined range.


EOB – MR LAVA

Liver intensity is represented by the last peak of the histogram, and its position can change among cases.



Liver segmentation Model-based approach



Reliable ROI for vein segmentation

7

Model registration – ITK code snippet

```
typedef itk::MeanSquaresImageToImageMetric < ImageType, ImageType > MetricType;
typedef itk::TranslationTransform < double, 3 > TranslateType;
typedef itk::NearestNeighborInterpolateImageFunction< ImageType, double > InterpolatorType;
typedef itk::RegularStepGradientDescentOptimizer OptimizerType;
typedef itk::ImageRegistrationMethod < ImageType, ImageType > RegistrationType;

metric = MetricType::New(); // Set metric
trans = TranslateType::New(); // Define registration transform
interpolator = InterpolatorType::New(); // Define interpolator
optimizer = OptimizerType::New(); // Define optimizer
registration = RegistrationType::New(); // Define registration method

registration->SetMetric( metric ); // Set metric
registration->SetTransform( trans ); // Set transform
registration->SetInitialTransformParameters( trans->GetParameters() );
registration->SetInterpolator( interpolator ); // Set interpolator
registration->SetOptimizer( optimizer ); // Set optimizer
registration->SetFixedImage( fix_image ); // Set fixed image (clustered input image)
registration->SetMovingImage( mov_image ); // Set moving image (liver model)
registration->SetFixedImageRegion( fix_image->GetLargestPossibleRegion() );

optimizer->SetMaximumStepLength( 10.0 ); // Set optimizer parameters
optimizer->SetMinimumStepLength( 1 );
optimizer->SetNumberOfIterations( 25 );

registration->Update(); // Start the registration
// Final transform = registration->GetTransform()
```

8

Liver segmentation – ITK code snippet

```
typedef itk::FastMarchingImageFilter< ImageType, ImageType > FastMarchingFilterType;
typedef FastMarchingFilterType::NodeContainer NodeContainer;

// Starting region is defined by the core of the registered model
...
NodeContainer::Pointer alives = NodeContainer::New(); // Internal voxels of starting region
NodeContainer::Pointer trials = NodeContainer::New(); // Contour voxels of starting region
...
// Speed of contour propagation is defined based on intensity similarity to starting region and
// probabilities in the registered model
...
ImageType::Pointer speed = ...;
FastMarchingFilterType::Pointer fastMarching = FastMarchingFilterType::New();


fastMarching->SetAlivePoints( alives ); // Set starting region
fastMarching->SetTrialPoints( trials );

fastMarching->SetInput( speed ); // Set speed function
fastMarching->SetStoppingValue( stop_value ); // Set maximal number of iterations


fastMarching->Update(); // Start the segmentation
// Result = fastMarching->GetOutput();
```

9


Liver segmentation Results – normal cases



Liver segmentation Results – cirrhotic cases



Liver segmentation Results – low contrast cases



Vessel segmentation

Vessel segmentation Hessian filtering

$$\frac{\partial}{\partial x} L(x, s) = s^7 L(x) * \frac{\partial}{\partial x} G(x, s)$$

$$|\lambda_1| \leq |\lambda_2| \leq |\lambda_3|$$

2D		3D		orientation pattern	
λ_1	λ_2	λ_1	λ_2	λ_3	
N	N	N	N	N	noisy, no preferred direction
		L	L	H-	plate-like structure (bright)
		L	L	H+	plate-like structure (dark)
L	H-	L	H-	H-	tubular structure (bright)
L	H+	L	H+	H+	tubular structure (dark)
H-	H-	H-	H-	H-	blob-like structure (bright)
H+	H+	H+	H+	H+	blob-like structure (dark)

Vessel segmentation Frangi's multi-scale vessel enhancement filtering

$$\mathcal{R}_A = \frac{(\text{Largest Cross Section Area})/\pi}{(\text{Largest Axis Semi-length})^2} = \frac{|\lambda_2|}{|\lambda_3|}$$

$$\mathcal{R}_B = \frac{\text{Volume}/(4\pi/3)}{(\text{Largest Cross Section Area}/\pi)^{3/2}} = \frac{|\lambda_1|}{\sqrt{|\lambda_2\lambda_3|}}$$

$$S = \|\mathcal{H}\|_F = \sqrt{\sum_{j \leq D} \lambda_j^2}$$

$$\mathcal{V}_o(s) = \begin{cases} 0 & \text{if } \lambda_2 > 0 \text{ or } \lambda_3 > 0, \\ (1 - \exp(-\frac{\mathcal{R}_A^2}{2\sigma^2})) \exp(-\frac{\mathcal{R}_B^2}{2\beta^2})(1 - \exp(-\frac{s^2}{2\sigma^2})) & \end{cases}$$

$$\mathcal{V}_o(\gamma) = \max_{s_{\min} \leq s \leq s_{\max}} \mathcal{V}_o(s, \gamma)$$

Vessel segmentation Based on vessel-enhancement

Vessel enhancement – ITK code snippet

```

// Create Hessian filter
typedef itk::HessianRecursiveGaussianImageFilter< ImageType > HessianFilterType;

HessianFilterType ::Pointer hessianFilter = HessianFilterType ::New();
hessianFilter->SetInput( input ); // Set input (MR image)
hessianFilter->SetSigma( sigma ); // Set radius of the tubular structure to enhance
hessianFilter->Update();

// Create vesselness filter based on "3D Multi-scale line filter for segmentation and visualization
of curvilinear structures in medical images", Yoshinobu Sato et al.
typedef itk::Hessian3DToVesselnessMeasureImageFilter<float> VesselnessFilterType;

VesselnessFilterType::Pointer vesselnessFilter = VesselnessFilterType ::New();
vesselnessFilter->SetInput( hessianFilter->GetOutput() );
vesselnessFilter->Update();

// Create intensity rescale filter to normalize vesselness values
typedef itk::RescaleIntensityImageFilter< FloatImageType, ImageType > RescaleFilterType;

RescaleFilterType ::Pointer rescale = RescaleFilterType::New();
rescale->SetInput( vesselnessFilter->GetOutput() );
rescale->SetOutputMaximum( 1000 );
rescale->SetOutputMinimum( 0 );
rescale->Update(); // Result is returned by rescale->GetOutput();
    
```

Vessel segmentation Results – Good contrast (40% of cases)

Good result

Vessel segmentation
Results – Low contrast (30% of cases)

Under-segmented result

19

Vessel segmentation
Results – Cirrhotic cases (30% of cases)

Under-segmented result

20

Challenges

21

Hessian vessel enhancement
Challenge 1: disconnected vessel tree

Vesselness can be nearly zero at junctions

22

Hessian-based vessel enhancement
Idea: Incorporate vessel orientation

23

Segmentation approach
Fuzzy connectedness*

- Theoretical background
- Graph-based approach based on the similarity of neighboring voxels
- Computes the connected components finding the optimal routes
- Advantages
 - Extend concept: incorporate local direction vector when similarity is computed (+ intensity and vesselness)
 - Fast and easy to implement


*Fuzzy Connectedness and Image Segmentation, J. K. Udupa, P. K. Saha, Proceedings of the IEEE, VOL. 91, NO. 10, OCTOBER 2003

24

Proposed method

Combine intensity with vessel direction


- Pre-process
 - Compute vesseness and local vessel directions
 - Identify the most likely vessel voxels (starting points)
- Iterative segmentation
 - Incorporate intensity difference of neighboring voxels
 - Allow **faster contour propagation in local vessel direction** (slow in perpendicular directions)
 - Propagate local vessel direction** from starting points into less confident voxels of its neighborhood
- Post-process
 - Morphological filtering and connected components



25

Compute vessel orientation

Preliminary results




Vectors located at 2 sides of a disjoint branch point in the same direction

26

Segmentation based on intensity and directions

Preliminary results




Promising results using one scale only

27

Multi-scale vessel enhancement

Challenge 2: leak between branches




Larger scales blur the small vessels

28

Scale selection

Preliminary results I.




Difficult to combine the information of different scales

29

Scale selection

Preliminary results II.





Difficult to combine the information of different scales

30

Future work and open questions


- Method
 - How to separate confident directions from noise?
 - How to combine the information of different scales?
 - Vector propagation to all neighbors or only in local direction?
 - What is the optimal balance between intensity and direction?
- Validation
 - Quantitative evaluation on synthetic images
 - Qualitative evaluation on real images
- Optimization
 - Reduce running time (1.5m >> 10-15s)
 - Reduce extreme memory usage (24 x input size)

Preliminary results


31

Improvement opportunities




Lesion detection and segmentation


Liver segmentation and morphological closing




Intensity normalization based on liver mode




Detect and segment lesions



Remove lesions from vessel









33

Normalize intensity heterogeneity

Normalize intensity based on the difference of the local average and the liver mode












34


Adaption to CT modality











35

